



OHAUS Electronic Balances

- Portable
- Scout Series
- Precision

- Affordable, superior performance from the market leader
- Large, high contrast LCD display for great visibility
- Quick and easy to use with large weighing pan and simple two button operation
- Customized performance with multiple weighing units, parts counting and battery operation

Standard Features: Large, high contrast LCD display, large stainless steel weighing surface, low profile, ergonomic design, integral security bracket, auto shut-off (user enabled/disabled), weigh below hook, sealed front panel and molded spill ring, stability indicator, error indicator for over/under loads, low battery indicator, span calibration mass.

Optional Accessories: Security device, carrying case, linearity calibration masses, scoops.

Industry Leading Warranty: All Ohaus Scout Balances are protected by a 5 year warranty on labor and parts. Our rugged construction and stringent quality control have been hallmarks of Ohaus scales and balances since 1907.



Ordering Information

ITEM	DESCRIPTION
SC2020	200g x 0.01g
SC4010	400g X 0.1g
SC6010	600g X 0.1g

Specifications

	SC2020	SC4010	SC6010
Capacity (g)	200	400	600
Readability (g)	0.01	0.1	0.1
Weighing modes	g, oz, oz t, dwt, parts counting		
Repeatability (Std.dev.) (g)	0.01	0.1	0.1
Linearity (g)	±0.01	±0.1	±0.1
Tare range	To capacity by subtraction		
Over range capacity	103%		
Stabilization time (seconds)	3		
Operating temperature range	50° to 104°F/10° to 40°C		
Calibration	Auto calibration		
Power requirements	100, 120, 220, 240 V ac, 50/60 Hz Adaptor (supplied) or one 9V battery (not included)		
Display (in/cm)	LCD (0.7/1.8 high)		
Pan size (in/cm)	4/10.2 dia	5x5.75/	5x5.75/
		12.7x14.6	12.7x14.6
Dimensions	6.75x2.4x8.25/17.1x6.1x21		
Net weight (lb/kg)	1.5/0.7	1.8/0.8	1.8/0.8
Shipping weight (lb/kg)	3.2/1.5	3.5/1.6	3.7/1.7