



## pH/Conductivity Meters

### pH/CON 10



- Toggle between pH/oC and conductivity/oC readings with one button press
- Includes multi-sensor submersible probe with 10-ft submersible cable

### pH/CON 300



- Waterproof and dustproof with IP67 rating-It FLOATS!
- Toggle between pH, TDS, and conductivity readings
- Dual display shows measured parameter plus temperature

#### Specifications

	pH/CON 10
<b>Range</b>	
pH	0.00 to 14.00 pH
Temp	0.0 to 100.0°C
Conductivity	0 to 19.99µS, 0 to 199.9µS; 0 to 1999µS, 0 to 19.99mS
<b>Resolution</b>	
pH	0.01 pH
Temp	0.1°C
Conductivity	0.01µS, 0.1µS, 1µS, 0.01mS
<b>Accuracy</b>	
pH	±0.01 pH
Temp	±0.5°C
Conductivity	±1% full scale or ±1 digit

#### Ordering Information

ITEM	DESCRIPTION
WD-35630-00	Standard pH/CON 10 meter
WD-35630-02	Waterproof pH/CON 10 meter
WD-35630-60	Standard pH/CON 10 meter kit
WD-35630-62	Waterproof pH/CON 10 meter kit
WD-35630-69	Kit without meter
WD-35630-50	Replacement pH/cond/°C probe with 10ft submersible cable
WD-35615-75	Belt loop portable meter carrying case.
WD-35615-07	AC adapter. Use with standard pH/CON 10 only.

#### Specifications

	pH/CON 300
<b>Range</b>	
pH	-2.00 to 16.00 pH
Conductivity	0 to 19.99µS, 0 to 199.9µS, 0 to 1999µS, 0 to 19.99mS, 0 to 199.9mS
TDS	0.00 to 9.99ppm, 10.0 to 99.9ppm, 100 to 9.99ppm
Temp	1.00 to 9.99 ppt, 10.0 to 9.99ppt, 100 to 200 ppt
Temp	0.0 to 100.0°C, 32.0 to 212.0°F
<b>Resolution</b>	
pH	0.01 pH
Conductivity	0.01µS, 0.1µS, 1µS, 0.01mS, 0.1mS
TDS	0.01, 0.1, 1 ppm; 0.01, 0.1, 1 ppt
Temp	0.1°C, 0.1°F
<b>Accuracy</b>	
pH	±0.01 pH
Conductivity	±1% full scale
TDS	±1% full scale
Temp	±0.5°C, ±0.5°F

#### Ordering Information

ITEM	DESCRIPTION
WD-35631-00	pH/CON 300 meter
WD-35631-60	pH/CON 300 meter kit
WD-35641-51	Replacement pH electrode
WD-35608-50	Replacement combination conductivity/TDS temperature probe
WD-35618-05	Temperature probe