

Portable Refractometer

Master M Series

Item Number: ATA-2313, ATA2323

- Easy user friendly portable digital wine refractometer
- Unique design for testing variety of samples
- Wide measurement range
- Easy to clean



The Mater M series of refractometers are new hand-held refractometers developed with high quality and ease of use in mind. The refraction view is brighter than that of all current models. The prism housing of the Master M refractometer is made of metal allowing for the temperature of hot or cool samples to conform to the temperature of the prism quickly. The Master M series is the successor model to the conventional N-1a .

Specification

Model	Master M	Master 2M
Part Number	ATA-2313	ATA-2323
Scale Range	Brix 0.00 to 32.0 %	Brix 28.0 to 62.0 %
Minimum Scale	Brix 0.2 %	Brix 0.2 %
Size and Weight	3.3 x 3.3 x 20.4cm, 160g	3.3 x 3.3 x 20.4cm, 160g